

TITLE

PROJECT

Continued from page

9/19 Repeat PCR dTAG 9/15 → Success

	1X
2x Phire	12.5
10 μ M Primer Mix	1
released DNA	1
dw	10.5
	<u>25</u>

9/20 Colony PCR PBA904 SOX9 g1, g2, g3, g4, NTCR3

	1X	25X
2x Phire MM	12.5	312.5
10 μ M each primer	1	25
Bacteria	1	25 -
dw	11	275
	<u>25</u>	

9/21 SOX9-g1 2nd RNP for DTAG:

1) mCherry Cassette	1	} BFP+
2)	2	
3)	3	
4) BFP Cassette	1	} mCherry +
5)	2	
6)	3	

Step 1: Prepare guide in MM → distribute to 8-strip

501281934	crRNA (100 μ M) SOX9-g1	1X	6X
1073191	tracrRNA (100 μ M)	2	12
		2	12

Incubate 95°C for 5 min → Cool to 4°C

Step 2: Prepare Cas9 RNP in 1.5-mL Eppendorf

	1X	6X
Cas9	2	12
PBS	3	18

Add 5 μ L Cas9 into 4 μ L guide RNA per rxn

Incubate Room Temp for 30 min

Continued to page

READ AND UNDERSTOOD BY:

SIGNATURE

DATE

SIGNATURE

DATE

TITLE

PROJECT

Continued from page

Step 3 Homology Cassette (read purified PCR product)
1000 ng needed

5	BFP mCherry Cassette (113 ng/μL)	1X	3X) add to RNP mix divided into 2 tubes
		8.84	26.5	
	BFP cassette (149 ng/μL)			
		6.71	20.1	

10	Lonza Master Mix	1X	6.5X
	P3 Solution	16.4 μL	106.6
	Supernatant PSI solution	3.6 μL	23.4

Step 4 : Prepare Cells

- 15
- Add 2 mL 0.25% Trypsin-EDTA to 10-cm T129 (mCherry + / BFP+)
 - Quench in 6 mL media
 - Conc:
 - Spin down 2 min 400g _____ mL of cell suspension. Decant
- 20
- ↳ need 0.5 mil cells / sample
 - 3 samples = 1.5 million cells BFP+
 - 3 samples = 1.5 million cells mCherry

- 25
- Resuspend in Lonza Master Mix (20 μL / sample)
 - Distribute into 8-strip containing RNP mix (3 each)
 - Transfer to electrode container

Step 5: Lonza

- 30
- FF 139 program for T129 cell line on X-Unit

Step 6: Plate in 12-well plate

- 35
- Add 100 μL media to well in 96-well of electrode
 - Transfer to 12-well
 - Top off to 200 μL per well

Continued to page

READ AND UNDERSTOOD BY:

SIGNATURE

DATE

SIGNATURE

DATE